



Education – Information Checklist

- Where to place the PAP in the bedroom
- How to adjust and position the head gear and mask
- How to turn on the PAP unit
- How to activate the ramp function
- How to attach the mask, pillows, headgear and tubing
- How to change or clean filters
- How to clean PAP unit and accessories
- How to disinfect the humidifier chamber
- Understand follow up program
- Refer to the manufacturers manual provided with your machine for further information

Safety

- Importance of grounding all electrical equipment
- How to obtain help if a medical emergency arises
- How to reach us for questions

Cleaning Instructions

- For mask, tubing, and water chamber: clean twice per week using Ivory/Dove dish soap and water. Allow to completely dry prior to use again
- Water chamber cleaning- soak in 4-5 tablespoons white vinegar and warm water for 30 minutes then rinse
- Filter- flick the dust off of this at least weekly

Questions

Call us anytime at (330) 572-1011